



Enamel Trays

Make: LabTek
Model: SUN-006
Origin: India

We are an authentic supplier and trader of Enamel Tray. This product has wide usage in different laboratories in different Sizes as :-

Enamel Tray/Sample Tray-Sizes(inch)

6 x 8, 8 x 10, 10 x 12
12 x 14, 12 x 18 15 x 18
18 x 24, 22 x 26 22 x 32, 24 x 36

Sampling Bowls

Make: LabTek
Model: SUN-007
Origin: India
Specification:

Different sizes made up of stainless steel.

Sizes: 25 cm dia., 30 cm dia. and 40 cm dia.



Stands, Bossheads And Clamps

Make: LabTek
Model: SUN-008
Origin: India

Laboratory clamps consist of a shaft and bosshead secured with thumbscrews onto a supporting rod or scaffold, and adjustable closing jaws for gripping objects safely and securely.

Most are metal clamps for the laboratory, with chemical and heat resistance, although the adjusting screws can be plastic, and the jaws may be lined with cork, plastic or rubber for enhanced grip.

- 🍏 General purpose plastic clamps
- 🍏 Retort clamps for round bottom-flasks
- 🍏 Glass joint clamps for ground glassware
- 🍏 Three-prong extension clamps for irregular or outsized objects
- 🍏 Chain clamps, screen support clamps, and other specialised clamps for burettes, condensers or reaction vessels

Laboratory Scoops

Make: LabTek
Model: SUN-009
Origin: India

Capacity 325 cc, 500 cc, 1000 cc, 2000cc and 5000cc

Material of Construction Stainless Steel, Aluminium

Finishing Mirror polishing



Pipette And Burettes

Make: LabTek
Model: SUN-010
Origin: India

Model No.	Description	Capacity
SUN-010-A	Graduated Pipette	10 ml (nominal)
SUN-010-B	Graduated Pipette	20ml (nominal)
SUN-010-C	Bulb Pipette	50m
SUN-010-D	Burette	50ml (nominal)
SUN-010-E	Burette	100ml (nominal)

