

## Laboratory Thermometers

Make: LabTek

Model: SUN-029

Origin: India

General Purpose Liquid-in-Glass Partial Immersion Thermometers are ideal for use in chemical laboratories, universities, scientific institutions, food and beverage processing, petrochemical laboratories and the pharmaceutical industry.

🍏 Material Glass

🍏 For Use With (Application) Chemical laboratories, universities, scientific institutions, food and beverage processing, petrochemical laboratories, pharmaceutical industry, medical industry, manufacturers, dairies, breweries and distilleries

🍏 Immersion Type Partial

🍏 Thermometer Fill Blue Organic Liquid



## Room Thermometer

Make: LabTek

Model: SUN-031

Origin: India

Types of room thermometer are normal room thermometer.

Range -30 to +60°C



## Digital Temperature Humidity Meter

Make: LabTek

Model: SUN-033

Origin: India

Large character LCD display, easy to read  
Temperature, Humidity and Clock display simultaneously  
Hourly chime, Daily alarm function

Max/Min temperature and humidity memory function

🍏 Calendar display function (Month and date)

🍏 Temperature display modes: °C/°F

🍏 2 time display modes: 12/24H

🍏 Desktop place or Wall hanging

Product Size (HxWxD): 105 × 97 × 23 mm

## Dial Thermometer

Make: LabTek

Model: SUN-030

Origin: India

DIAL THERMOMETERS work on principle of bi metal expansion. Scale or range of temperature which DIAL THERMOMETER can measure may vary from -200 0 C to 650 0 Centigrade. DIAL THERMOMETER consists of a sharp metal indicator which is fixed on a circular temperature scale. There are two main components, one is a needle and another is a temperature sensor. Few more types of DIAL THERMOMETERS are Liquid type, Gas filled and Va pour tension dial thermometers.



## Maximum & Minimum Thermometers

Make: LabTek

Model: SUN-032

Origin: India

A maximum–minimum thermometer. The scales are Fahrenheit on the inside of the U and Celsius on the outside. The current temperature is 23 degrees Celsius, the maximum recorded is 25, and the minimum is 15, both read from the base of the small markers in each arm of the U tube.

🍏 Ranges Available: -30° C to +50 ° C



## Electronic Thermometer

Make: LabTek

Model: SUN-034

Origin: India

This series of thermocouple thermometers are micro controller based testers. With very intuitive interface, the user can fully operate the tester without second thought. They are digitally calibrated for long time stability. The whole series are equipped with PC interface for advance applications.

Features:

🍏 Timer function

🍏 PC interface

🍏 Auto power off

🍏 Resolution 0.1°C; 0.1°F

Backlit LCD