

Sampling Trays

Make: LabTek

Model: SL-001

Origin: India

These original white polypropylene trays are for field work and general lab use. Polypropylene is resistant to most chemicals and solvents and can be autoclaved.

The small trays are great for pond sampling and the larger ones for sorting leaf litter. We do not recommend the larger tray for aquatic sampling as it is a little too flexible to be carried full of water.



Measuring Tapes

Make: LabTek

Model: SL-002

Origin: India

Material: Plastic

Colour: Sky Blue & Black

Size (L x B x H): Tape: 10 cm x 4 cm x 9 cm; Level: 6 cm x 3 cm x 6 cm

No. of Pieces: 2

Included Components: 2 Measuring Tapes

Brand Name: FREEMANS

Colour: Sky Blue & Black, Sky Blue & Black

Included Components 2 measuring tapes

Item Weight 499 grams

Material Plastic


Graduated Measuring Jars

Make: LabTek

Model: SL-003

Origin: India

These are commonly used at site for measurement of Silt & Sand, to measure the quantities of liquids these are available in Ordinary Glass, Polypropylene Plastic, Borosilicate Glass, Soda Glass. Capacity available from 5ml to 2000ml.



Funnels and Volumetric Flasks

Make: LabTek

Model: SL-004

Origin: India

1 Volumetric Flask, Capacity 100ml

2 Volumetric Flask, Capacity 500ml

3 Volumetric Flask, Capacity 1000ml

4 Conical Flask, Capacity 500ml (nominal)

5 Conical Flask, Capacity 250ml (nominal)

6 Filter Flask, Capacity 1000ml

7 Extraction Flask, Capacity 250ml

8 Buchner Funnel

9 Glass Funnel

10 Funnel

11 Funnel

Moisture Tins

Make: LabTek

Model: SL-004

Origin: India

To determine moisture contents in soil and other material. These are made from aluminium or stainless steel.

Sizes available are :

(a) 25 mm dia

(b) 50 mm dia

(c) 75 mm dia

(d) 80 mm dia

(e) 100 mm dia.

