

## Water Bath (6 Holes)

Make: LabTek

Model: SUN-011

Origin: India

Our boiling water baths come in rectangular shape and engineered to meet most demanding needs of heated applications. We make 3 models that are available in 4 (four), 6 (six) and 12 (twelve) holes. Temperature range is Ambient +5°C to 100°C, which is controlled through thermostatic regulator. Temperature can be seen on LED display. These water baths are constructed with 304 stainless steel for inner chamber; exterior using powder coated GI sheet. External and internal surfaces feature corrosion resistance property and also easy to wipe and clean after use. These units also have load indicators which indicate running status of machines. Water drain nozzle is fitted at back side of the machines.



Model	BH4	BH6	BH12
No of holes	4	6	12
Lid	Concentric rings 75 mm dia.	Concentric rings 75 mm dia.	Concentric rings 75 mm dia.
Temperature	Ambient+5 to 100°C	Ambient+5 to 100°C	Ambient+5 to 100°C
Temperature controller	Thermostatic rotary knob	Thermostatic rotary knob	Thermostatic rotary knob
Display	LED Display	LED Display	LED Display
Construction	Double walled	Double Walled	Double walled
Inner MOC	SS 304	SS 304	SS 304
External MOC	Powder coated GI sheet	Powder coated GI sheet	Powder coated GI sheet
Insulation	Glass wool	Glass wool	Glass wool
Power supply	220 Volts 50 Hz	220 Volts 50 Hz	220 Volts 50 Hz
Optional	<ul style="list-style-type: none"> <li>Digital temperature controller</li> <li>SS 304 exterior</li> </ul>	<ul style="list-style-type: none"> <li>Digital temperature controller</li> <li>SS 304 exterior</li> </ul>	<ul style="list-style-type: none"> <li>Digital temperature controller</li> <li>SS 304 exterior</li> </ul>



## Universal Water Bath Serological

Make: LabTek

Model: SUN-012

Origin: India

Matching the various IS standards are Shaker Water Bath made by a classic that reciprocates shaking motion application in the ever growing field related to agriculture, textile, medical, industry, research laboratories and hospitals.

- 🍏 Material : Stainless Steel
- 🍏 Shape : Rectangular
- 🍏 Ambient : + 5 degree C to 99 degree C
- 🍏 Motor: AC/DC Gear Motor
- Accuracy: +/- 0.5 degree C
- Revolutions per minute 50 To 150

## Universal Water Bath with Stirrer

Make: LabTek

Model: SUN-013

Origin: India

Matching the various IS standards are Shaker Water Bath made by a classic that reciprocates shaking motion application in the ever growing field related to agriculture, textile, medical, industry, research laboratories and hospitals.

- 🍏 Material: Stainless Steel
- 🍏 Shape Rectangular
- 🍏 Ambient + 5 degree C to 99 degree C
- Motor AC/DC Gear Motor
- Accuracy +/- 0.5 degree C
- Revolutions per minute 50 To 150

