

Ultra Low Temperature Laboratory Deep Freezer MODEL - SL - 0970, MAKE - SLBT, ORIGIN - INDIA



Product Description

These Low Temperature Freezers are manufactured using sophisticated plant machinery that are coupled with advanced technology, to ensure high standards of quality. These products are offered with thermostat or digital controller for effective monitoring of temperatures of these products. The inner walls of these products are manufactured from high quality acrylic to avoid temperature loss.

Specifications:

Temperature ranges: -20°C to ambient / -40°C to ambient / -60°C to ambient / -80°C to -40°C

Controlling accuracy +2°C @ Minimum.

Display accuracy +0.5%, full scale.