

Vortex Mixers

MODEL - SL - 1036, MAKE - SLBT, ORIGIN - INDIA



Specification

Material	Mild Steel
Automation Type	Automatic
Application	Shaker
Motor	Brushless DC
Voltage	220 V
Wattage	15 Watts

This **Vortex Mixer** is vigorous agitating equipment used for thorough mixing of precipitates in test tubes, small bottles or flasks. The motor driven neoprene cup mounted eccentrically on the shaft vibrates rapidly. Varying the angle and pressure applied controls the mixing action. Provided switch to select direct working or on touch working. A speed regulator switch also provided to control the speed.