


Notch Apparatus

MODEL – SL - 1200, MAKE - SLBT, ORIGIN – INDIA



Notch Tank :

Notch Tank of suitable size having arrangement to fix replaceable notches with a set of welded mesh baffles for accommodating $\frac{3}{4}$ "1" crusher stones for steadying the flow of water Notches- (Brass): a. Trapezoidal Notch b. Rectangular Notch-100mm width

Hook Gauge:

30cm range with mounting provision

Measuring Tank:

M.S, Power Coating/M.S.FRP Lining/S.S Unit of auditable size with overflow arrangement and provided with a PVC control valve and an elbow.

Sump Tank:

M.S. Powder Coating /MS.FRP Lining/SS. Unit of suitable size for storage of water to circulate through experimental unit.

Supply Pump Set:

0.5 HP single phase [ISI marked] to pump water from the sump from the sump fitted with a control valve

Main switch :

Single phase main switch

Frame Work Mounting:

Mounted on frame work and completely fitted with all the above items as self contained unit suitable for operation without and foundation.

Services Required:

Single phase power supply 200/220V and water to fill up sump Tank.