

## Flame Photometer

**MODEL – SL - 1258, MAKE – SLBT, ORIGIN – INDIA**



Display	Graphical LCD 128x64
Portable	No
Application	Industrial Use
Automation Grade	Automatic
Capacity	0.5 ml.
Color	White
Product Type	Flame photo meter
Power Supply	Smps Power Supply
Filters	Na & K Supplied (Ca & Ba are optional at extra cost)
Atomizer	Axial Flow Type
Average Time	in-build Software
Calibration	upto 5 Point With curve fitting software
Curve Accuracy	<2% error when 3 ppm Na/K and 5 ppm, Li are set as maximum Standards
Detection	Silicon Photodiode
Flame System	LPG and Dry Oil Free Air
Warranty	12 Months
Gas Control	Adjustable with knobs
Ignition System	Auto Ignition System
Resolution	0.1
Sensitivity	Na: 0.5 ppm, K: 0.5 ppm

