

**Automatic Fat Analyzer**  
**MODEL – SL - 1492, MAKE – SLBT, ORIGIN – INDIA**



**Technical parameters**

- Measuring range: oil content within the scope of 0.5%-60% of good, feed, oil and various kinds of fat products
- Measuring sample quantity: 6 simultaneously
- Heating time: within 10 min
- Recovery system: automatic
- Recovery rate:  $\geq 80\%$
- Temperature controlling range: room temperature 0-100°C
- Power: AC220V/50Hz, 1000W