

Chemical Storage Cabinet Make SLBT Model SL-1549 Origin India



Specifications

Robust construction and longevity - body made of sheet steel with high-grade structured paint finish.

No unauthorized use - doors lockable with cylinder locking (integrated in the turning handle)

Easy alignment - adjusting aids to compensate for uneven floor

Reduction of odors inside the cabinet - natural ventilation due to vent openings at the cabinet bottom

Shelves - available either in galvanized or painted version, light grey, height adjustable, load capacity 100 kg (evenly distributed load)

Pullout shelves - robust, high-quality pullout rails with pullout shelves made of metal and removable plastic sump, load capacity 50 kg (evenly distributed load)