

Metric Weighing Set Make SLBT Model SL-1675 Origin India



Product Specification

This metric weight set consists of high quality cast and brass metal weights. The set has a total of 10 weights and includes the following pieces: one piece of 1kg one piece of 500g two of 200g one of 100g one of 50g two of 20g one of 10g one of 5g The total mass of the set is 2105 grams and it provides for weighing objects from 5g to 2105g with 5 grams increments. This set is perfect in combination with a balance scale, to accurately measure the mass of the object with an unknown mass.